

Laser Practice Eye

DEVICE DESCRIPTION

- The optics resembles an aphakic, emmetropic eye.
- A photographic film is used as an artificial diseased retina, which responds to the laser shots. Under exposure / correct exposure and over exposure gives a different visible reaction.
- It can be mounted on slit lamp of all available makes.
- Laser Practice Eye is an optical + mechanical device which does not require electricity.

INTENDED USE// INDICATION FOR USE

Laser Practice Eye is an ophthalmic teaching device used for simulating green/red laser shots in an artificial diseased retina.

CONTRAINDICATIONS

Not Known

HOW SUPPLIED

Laser Practice Eye is supplied un-sterile, in individual leather box. The overall packing contains Laser Practice Eye and medical literature insert, in an outer leather box, with external identification stickers.

⚠ WARNING

- The power setting used to get a reaction in the Laser Practice Eye may be higher than required in the human eye. The instructor needs to inform the trainees about this fact.
- Only qualified medical practitioner should use the device.
- Read instruction before use.
- Store at room temperature.
- For External use only.
- Keep out of reach of children.

PRECAUTIONS

- When handling the product care should be taken to avoid damage from handling.

DIRECTION FOR USE

- Mount the Laser Practice Eye on the chin rest of the Slit Lamp.
- Load the new photograph of diseased retinal image in the eye by unscrewing the back cover.
- Set the parameters of the Laser as per manufactures instructions.
- Place the desired contact laser lens with a coupling medium on the cornea.
- Focus the aiming beam through the contact laser lens before firing the laser beam.
- Adjust power of the laser beam till a visual reaction can be easily observed.
- Remove the photograph to view the reaction achieved.



DISPOSAL OF THE USED DEVICE

The used device should be disposed off in compliance with the local regulatory requirements.

REUSE OF DEVICE

The device should undergo a detailed multi step process of cleaning & sterilization before reuse.

RETURN OF DAMAGED PRODUCT

Return the product in its original packing identified by the batch number, purchase information, your reference and reason for return. Please contact your local distributor office regarding product return/exchange.

EXPLICATION OF SYMBOL:

Symbol	Meaning
	Caution
	Manufacturer
	Date of Manufacture
	Keep away from sunlight
	Consult instructions for use
	Re-use
	Authorized Representative in the European Community
	This side up
	Do not use hook

Symbol	Meaning
	Conformity of European Norm
	Batch Code
	Temperature Limit
	Do not use if package is damaged
	Humidity Limitation
	Handle with care
	Maximum stacking
	Medical Device
	Non-Sterile

MADHU INSTRUMENTS PVT. LTD.

Mfg. Facility: A 260, Okhla Industrial Area Phase-I, New Delhi-110020 India.
 Regd. Add.: F-90/3D, Okhla Industrial Area Phase-I, New Delhi-110020 India
 Customer Care No.: +91-11-42701030, 42701031
 Email: customercare@madhuinstruments.com
 Website: www.madhuinstruments.com

EC REP

Obelis s.a.
 Bd. Général Wahis 53
 B-1030 Brussels, Belgium
 Phone: 32.2.732.59.54, Fax: 32.2.732.60.03
 E-mail: mail@obelis.net