



## Capsular Support Segment

This is a partial ring of clear PMMA covering approximately one quadrant with a scleral fixation eyelet for temporary or permanent fixation.



### Product Description

Product Code	Model	Design	Normal Size	Compressed Size	Eyelets	Material	Qty. per box
MIPL/D11	CS Segment (Type 10)		–	–	3	PMMA	1 pc.

### Advantages

- It is used to stabilize and recenter the capsular bag
- Specially designed for scleral fixation with sutures through eyelets in a cataract surgery.
- The eyelet is projecting anteriorly which allows it to sit outside the capsular bag
- CE certified

## Indications

- Cases with extreme zonular weakness
- Cases with missing zonules larger than 4 clock hours
- Lens subluxation

## Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate . The international symbols used on the packaging are also explained.

## Sterilization

Products are sterilized through Ethylene Oxide (EO) according to a validated cycle.



**MADHU INSTRUMENTS PVT LTD**

A-260, Okhla Industrial Area, Phase-1, New Delhi -110020, India

T: (+91-11) 4270 1030, 31 | E: Sales@madhuinstruments.com | W: www.madhuinstruments.com

Customer Care No. +91 9810866176 