



Backflush Handle (Active Aspiration)

This instrument incorporates a silicon backflush reservoir on the conventional Charles flute handle which helps in safe backflushing of the incarcerated tissue during active aspiration of intraocular fluids.



Product Description

Product Code	Description	Qty. per box
MIPL/E17	Backflush Handle (Active Aspiration)	1 pc. (Sterile)

Features

- EO Sterile
- Disposable

Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate . The international symbols used on the packaging are also explained.

Sterilization

Products are sterilized through Ethylene Oxide (EO) according to a validated cycle.



MADHU INSTRUMENTS PVT LTD

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)