


Capsular Tension Ring with Scleral Fixation Arm

This is a modified CTR which has an additional fixation eyelet that allows it to be sutured to the scleral wall without piercing the capsular bag. This ring is used in cases where a standard CTR may not provide adequate centration and stability.



Product Description

Product Code	Model	Design	Normal Size	Compressed Size	Eyelets	Material	Qty. per box
MIPL/D3/T6	CTR Type 6 (L)		12.0 mm	10.0 mm	3	PMMA	1 pc. (Sterile)
MIPL/D3/T7	CTR Type 7 (R)		12.0 mm	10.0 mm	3	PMMA	1 pc. (Sterile)
MIPL/D3/T8	CTR Type 8 (LR)		12.0 mm	10.0 mm	4	PMMA	1 pc. (Sterile)

Advantages

- It is used to stabilize and recenter the capsular bag
- Specially designed for scleral fixation with sutures through eyelets in a cataract surgery.
- The eyelet is projecting anteriorly which allows it to sit outside the capsular bag
- A double eyelet design is also available for cases which require additional support
- CE certified

Indications

- Cases with extreme zonular weakness
- Cases with missing zonules larger than 4 clock hours
- Lens subluxation

Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate . The international symbols used on the packaging are also explained.

Sterilization

Products are sterilized through Ethylene Oxide (EO) according to a validated cycle.



MADHU INSTRUMENTS PVT LTD

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)