

Cellulose Sponge Spears

High quality cellulose sponge spears specifically designed for fast and high capacity fluid absorption.



Product Description

Product Code	Description	Qty. per box
MIPL/B1/5	Cellulose Spears per pouch, 50 pouches per box	250 Spears per box (Sterile)

Features

- Special cellulose material gives the fastest wicking spears available
- Highly absorbent and an ideal replacement for cotton swabs
- CE certified



Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate. The international symbols used on the packaging are also explained.

Sterilization

Products are sterilized through gamma radiation according to a validated cycle.



MADHU INSTRUMENTS PVT LTD

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)