


Cortex Aspirating Cannula

Disposable ophthalmic cannula used for aspirating cortical debris when chamber is maintained with either an anterior chamber maintainer or a viscoelastic substance.



Product Description

Product Code	Description with Gauge	Needle Length	Qty. per box
MIPL/E10/A1	Cortex Aspirating Cannula (Simcoe) - 21G	22 mm	10 Pcs. (Sterile)
MIPL/E10/A2	Cortex Aspirating Cannula (Simcoe) - 23G	22 mm	10 Pcs. (Sterile)
MIPL/E10/A4	Cortex Aspirating Cannula (Simcoe) - 25G	22 mm	10 Pcs. (Sterile)
MIPL/E10/A5	Cortex Aspirating Cannula (Simcoe) - 26G	22 mm	10 Pcs. (Sterile)
MIPL/E10/A3	Cortex Aspirating Cannula (J-Shaped) - 25G	22 mm	10 Pcs. (Sterile)



Features

- Thin wall cannula used for maximum flow
- CE certified

Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate. The international symbols used on the packaging are also explained.

Sterilization

Products are sterilized through Ethylene Oxide (EO) according to a validated cycle.



MADHU INSTRUMENTS PVT LTD

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)