

I – Stand (Eye Ball Stand)

This mechanical device is used for securely fixing the human or animal eyes and micro surgical training on the same.



Product Description

Product Code	Description	Qty. per box
MIPL/M1	Eyeball Stand	1 pc.

Features

- Both human and animal eyes can be fixed on the stand
- Required anterior chamber pressure can be created and maintained
- Can be sterilized by autoclaving or EO gas
- CE Certified

Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate. The international symbols used on the packaging are also explained.

Cleaning & Disinfection

For cleaning and disinfection of the device, kindly refer to the IFU supplied with the device or contact our customer care number



MADHU INSTRUMENTS PVT LTD

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)