


## Iris Retractors

Iris retractors are used for temporary retraction (pupil enlargement) of the iris where dilation cannot be achieved pharmacologically.



## Product Description

| Product Code | Description     | Qty. per box     |
|--------------|-----------------|------------------|
| MIPL/D1      | Iris Retractors | 5 Pcs. (Sterile) |

## Features

- Flexible retractors made of bright blue PP (Polypropylene)
- Stretches and stabilizes the Iris during the cataract surgery
- An adjustable silicon stopper provides a superior grip and allows adjustment of the pupil size
- Easy to insert and remove
- CE certified

## Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate. The international symbols used on the packaging are also explained.

## Sterilization

Products are sterilized through Ethylene Oxide (EO) according to a validated cycle.



**MADHU INSTRUMENTS PVT LTD**

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)