

Laser Practice Eye

Laser-i is an ophthalmic teaching device used for simulating green/red laser shots in an artificial diseased retina.



Product Description

Product Code	Description	Qty. per box
MIPL/M3	Laser Practice Eye	1 Pc.

Features

- The optics resemble an aphakic, emetropic eye
- A photographic film is used as an artificial diseased retina, which responds to the laser shots
- It can be mounted on slit lamp of all available makes
- "Laser-i" is an optical + mechanical device which does not require electricity
- CE Certified

Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate. The international symbols used on the packaging are also explained.

Cleaning & Disinfection

For cleaning and disinfection of the device, kindly refer to the IFU supplied with the device or contact our customer care number



MADHU INSTRUMENTS PVT LTD

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)