

## Vitreotomy Lens Set

This vitrectomy lens set comprises a set of 7 contact lenses and 2 rings to hold the contact lenses in position. The main advantage of this system is that it minimizes the dependence of a well-trained assistant.



### Product Description

Lens	Central Fundus	Biconcave (-90D)	Prism 20°	Prism 30°	Biconcave Prism	Magnifying Lens	Widefield
Code No.	MIL 1	MIL 2	MIL 3	MIL 4	MIL 4A	MIL 5	MIL 5A
Cross Section View							

## Features

- High quality glass optics for clear images
- Unique silicon storage cum sterilization container
- Autoclavable Lenses
- Suturable metal ring
- Silicon ring for short procedures
- Flexibility to order individual lenses and rings

## Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate . The international symbols used on the packaging are also explained.

## Sterilization

Products are sterilized through Ethylene Oxide (EO) according to a validated cycle.



**MADHU INSTRUMENTS PVT LTD**

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)