

## 90D Aspheric Lens

This slit lamp lens is used for general examination of the fundus and small pupil examinations.



### Product Description

Product Code	Description	Field of View	Qty. per box
MIPL/K3	90D Aspheric Lens	94°	1 pc.

### Features

- This Slit lamp lens is the most popular power for non-contact fundus examination.
- It offers a wide-field fundus imaging out to and beyond the equator.
- A high magnification provides an increased illuminated field area which makes it ideal for posterior pole scanning.
- Made in glass.
- Aspheric Optics
- CE Certified.

## Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate . The international symbols used on the packaging are also explained.

## Sterilization

Products are sterilized through Ethylene Oxide (EO) according to a validated cycle.



**MADHU INSTRUMENTS PVT LTD**

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)