


## Four Mirror Goniolens

This is a standard gonio lens which is used for static and dynamic gonioscopy.



### Product Description

Product Code	Description	Field of View	Qty. per box
MIPL/I4	Four Mirror Goniolens	62°	1 pc.

### Features

- Four mirrors angled at 62° help in viewing all quadrants of the anterior chamber without rotation.
- The small contact surface of 9 mm diameter makes the lens very convenient to use.
- It has a unique detachable handle.
- Coupling fluid not required.
- Made in glass.
- No Contact medium required.
- CE Certified.

### Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate. The international symbols used on the packaging are also explained.

### Cleaning & Disinfection

For cleaning and disinfection of the lens, kindly refer to the IFU supplied with the lens or contact our customer care.



**MADHU INSTRUMENTS PVT LTD**

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)