

## Two Mirror Goniolens

This gonio lens is used for viewing the anterior chamber angle and central retina.



### Product Description

Product Code	Description	Field of View	Qty. per box
MIPL/I2	Two Mirror Goniolens	62°/62°	1 pc.

### Features

- Two opposing mirrors inclined at 62° provide a complete 360° view of angle with only 90° of rotation of the lens.
- The small contact surface helps in examination of children and adults with small fissures.
- Made in glass
- Contact medium required
- CE Certified

### Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate. The international symbols used on the packaging are also explained.

### Cleaning & Disinfection

For cleaning and disinfection of the lens, kindly refer to the IFU supplied with the lens or contact our customer care.



**MADHU INSTRUMENTS PVT LTD**

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)